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Luminol chemiluminescence reagent kit

Product Description

This kit can be used for HRP-labeled Western Blot and ELISA experiments. X-ray film or luminometer can be used to detect the light intensity or image. The light intensity of the labeled antigen will be qualitatively or quantitatively detected. This kit has revolutionized the traditional ECL substrate system, resulting in excellent detection limits, duration of luminescence and reagent stability. Overall performance evaluation parameters would be enhanced more than hundred times. This kit is as same as the world-renowned brand's product on sensitivity systems in actual performance. Store at 2-8°C

Name	Catlog Number	Solution A ml	Solution B ml	L O D	Same as the flowing brands		
					Brand	Name	Cat
LumiA	SB-WB001S	50	50	10 pg	Super Signal®	Pico	34080
	SB-WB001L	250	250				
ECL Plus	SB-WB003S	50	50	5 pg	Super Signal®	Picoplus	34580
	SB-WB003L	250	250				
Super	SB-WB002S	50	50	1 pg	Super Signal®	Dura	34095
	SB-WB002L	250	250				
LumiQ	SB-WB012S	50	50	0.1 pg	Merck Millipore	Immobilon®	WBKLS0500
	SB-WB012L	250	250		Bio-Rad	Clarity®	1705061
Lumi Best	SB-WB011S	50	50	fg	Super Signal®	Femto	34095
	SB-WB011L	250	250				

Materials needed

1. Secondary antibody incubated Western blot membrane.
2. TBST or PBST buffer.
3. Chemiluminescence detection equipment or X-ray film.
4. If using X-ray film you need cassette and plastic film.

How to use

1. Western Blot experiments were performed as usual and the secondary antibody incubation was completed. The membrane was washed three

times for 5 min each with TBST or PBST buffer.

2. Take a certain volume of liquid A and liquid B, the 1:1 mixing reserve (For example, a 6CM * 9CM PVDF membrane needs 400ul + 400ul liquid A + B, skilled experimenter can also add less volume).

3. Remove the TBST or PBST buffer in the container. The mixed liquid A + B were evenly leached on the membrane.

Troubleshooting

Problem	Cause	Solution
Yellow band		
Light-emitting time getting shorter	Secondary antibody concentration is too high	Reduce the secondary antibody concentration
Black band		
Weak signal or no signal	Too much HRP in the system	Further dilute the HRP-conjugate
	No target protein in the sample	Check up sample preparation
	Incomplete transferring	Extend the transfer time
	Second antibody degeneration	Make new secondary antibody solution
Background is too dark or appear non-specificity bands	Poor primary antibody specificity	Adjusted primary antibody concentration
	Incomplete blocking	Extend the blocking time
	Blocking solution degeneration	Make new blocking solution
	Overexposure	Reduce exposure time

Notice: Basic (catlog: sb-wb001), ECL Plus (catlog: sb-wb003), and Super (catlog: sb-wb002) have same protocol, but they have different sensitive: Super is the highest, basic is the lowest, ECL plus is in the middle. If you want to buy again, Please see the catlog number on the bottle or box!

Shanghai Share-bio Technology Co, Ltd

Address: Room 3321, No. 1200 Cangyuan Road, Minhang District, Shanghai

<http://www.share-bio.com>

Tel: 021-60533564