Western Blot Fast Stripping Buffer (CAT: SB-WB007)

This product ship in RT, and store in 2-8 degree, please use within 12 months.

Product feature

Save your time, effort, and sample----multi target protein can be test in just one blot, more than 10 target protein can be probed in just 1 blot. Reduced the labor and sample that needed in multi sample loading and blotting, also can reduce the artificial caused by multi loading.

Ready and easy to use---- only 10min is needed to finish a stripping, no need for extra equipment, no need to heat, probe can be done after a simple wash.


Mild formula, environment and user friendly----mild formula, suitable for PVDF and NC membrane, suitable for fluorescent, HRP and NIR technology. Less effect on protein sample, to ensure multi time detection of interest protein. Biology safe and user-friendly formula, can be used in open lab space.

Procedures

1. Remove the redundant reagent by PBST/TBST/PBS/TBS/or water wash;
2. Add enough stripping buffer to cover whole membrane. Shake on a rotate for 10-20min in room temperature. Stripping time can be adjust depending on the stripping efficiency on target band. 20 min is enough for most hard-to-remove antibody, for example Tubulin.
3. Discard stripping buffer, and wash with wash buffer for 5min. the efficiency can be checked at this step by re-develop in a long developing time.
4. Block with your block buffer, and ready for next probe.
5. we have checked the room temperature stripping time can be as long as 6 hours, no significant signal loss for weak signal protein band.

Cautions to user

1. Repeat use is not suggested since the pH value is critical for efficient strip, reuse will decrease the efficiency.
2. Both NC and PVDF membrane can be used, but we recommend PVDF membrane since PVDF membrane have a better protein reservation function.
3. Buffer can be heated to 50 degree is not satisfied with the stripping effect, and a longer stripping time can help.
4. Reduce your developing time if HRP is used in developing, a strong developing will leave a yellow mark on your membrane, and its irreversible.
5. PPE is required for your safety. **For research use only, not for use in diagnostic procedures.**
6. Low temperature may cause precipitation, placed at room temperature for half an hour or dissolved in 50 °C water before use.

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