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## All-in-One First-Strand Synthesis MasterMix (with dsDNase)

REF:SB-RT001 20 μl\*50 T;20 μl\*100 T

**Storage Condition: -20°C**

### Components

Component	Specification (50T)	Specification (100T)
All-in-One First-Strand Synthesis Master Mix	200 μl	400 μl
dsDNase	50 μl	2×50 μl
10× dsDNase Buffer	200 μl	200 μl
Nuclease-Free Water	1 ml	2×1 ml

### Description

All-in-One First-Strand Synthesis MasterMix (with dsDNase) is a reverse transcription system developed for first-strand cDNA synthesis with added convenience. It contains all the components necessary for first-strand cDNA synthesis, e.g. MonScript™ RTase III, reaction buffer, RNase Inhibitor, dNTPs, Oligo(dT)<sub>20</sub>VN and random primers, and requires only the addition of RNA template and water. MonScript™ RTase III is the 3rd generation reverse transcriptase developed from the M-MLV RTase. It possesses no RNase H activity and has elevated reaction temperature of up to 55°C , which significantly increased reverse transcription efficiency and resistance to complicated templates, e.g. high GC contents and massive secondary structures. The enzyme provides better specificity and yield and higher chance to obtain full-length cDNA. By using this premix, cDNA of up to 12 kb can be obtained within 15 min that is recommended to be used in downstream qPCR experiments.

High-efficiency dsDNase is used in this premix. This enzyme is heat-labile and can be heat-inactivated rapidly and irreversibly with high temperature. Thus removal of gDNA and reverse transcription reactions are finished in one tube with one experimental set-up. Different from the usually used DNase I, using dsDNase requires no addition of EDTA. This avoids possible damage to the RNA template, reduces the chance of contamination, and saves the experiment time.

### Method

#### First-strand cDNA synthesis

- Set up the following reaction mix on ice:

Reagent	Volume
Template RNAa	50 ng~1 μg
All-in-One First-Strand Synthesis MasterMix	4 μl
dsDNase	1 μl
10× dsDNase Buffer	2 μl
Nuclease-Free Water	To 20 μl

- High-quality RNA templates extracted by a kit is recommended.
- Mix gently and spin briefly.
- Incubate under 37°C for 2 min to remove the gDNA contamination.
- Incubate under 55°C for 15 min.
- Terminate the reaction by incubate under 85°C for 5 min.
- Put the cDNA on ice for downstream experiments or store at -20°C immediately.

### Note

The premix contains Oligo(dT)<sub>20</sub>VN and random primers, and is suitable for eukaryotic mRNA with Poly(A), prokaryotic RNA and eukaryotic rRNA/tRNA without Poly(A), but is not suitable for small RNA templates like miRNA.

### Related Products

REF	Product Name	Specification
SB-MR009	Trizol 总 RNA 提取试剂盒	100 ml
SB-Q204	Universal qPCR Master Mix	5*1 ml