

**Product Datasheet****V5 tag Antibody**

Catalog No: SB-AB0052

Reactivity: ALL

Isotype: Rabbit IgG


Applications: WB

**Information****UniProt ID:****All Names:** V5 tag**Form:** Liquid**Storage instructions:** Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.**Storage buffer:** pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.**Purity:** Affinity-chromatography**Immunogen:** Synthesized peptide**Molecular Wt.:****Application**

WB: 1:3000~1:10000

**Background**

The V5 epitope tag is derived from a small epitope (Pk) present on the P and V proteins of the paramyxovirus of simian virus 5 (SV5). The V5 tag is usually used with all 14 amino acids (GKPIPNPLLGLDST), although it has also been used with a shorter 9 amino acid sequence (IPNPLLGLD).



Western blot analysis of extracts from v5-tag fusion protein, using V5 tag antibody.

**For Research Use Only. Not For Use In Diagnostic Procedures.**