

## Product Datasheet

# MBP Antibody

Catalog No: SB-AB0029

Reactivity: All

Isotype: Rabbit IgG

Applications: WB

### Information

**UniProt ID:** P0AEX9

**All Names:** MalE; malJ; Maltodextrin binding protein; Maltose ABC transporter periplasmic protein; Maltose binding periplasmic protein; MMBP; Periplasmic maltose binding protein;

**Form:** Liquid

**Storage instructions:** Store at +4° C short term. Store at -20° C long term. Avoid freeze / thaw cycle.

**Storage buffer:** pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Purity:** Affinity-chromatography

**Immunogen:** A synthesized peptide


**Molecular Wt.:** 37 kDa

### Application

WB 1:1000~1:5000

### Background

MBP-Tag antibody is part of the Tag series of antibodies, the best quality in the research.



Western blot analysis of MBP expression in E.coli lysates.

**For Research Use Only. Not For Use In Diagnostic Procedures.**